

da Vinci Xi[®]
SURGICAL SYSTEM
da Vinci X[®]
SURGICAL SYSTEM

EndoWrist

FORCE BIPOLAR

Addendum to the Reprocessing Instructions Appendices

Supplemental instructions for the cleaning,
disinfection, and sterilization information
for reusable instruments and accessories used
with *da Vinci Xi*[®] and *da Vinci X*[®] Systems

The reprocessing instructions and parameters given represent
the manufacturer's recommendations according to
EN ISO 17664:2004 and ANSI/AAMI ST81:2004/(R) 2010.

This Addendum is intended to be used together with the
Reprocessing Instructions Appendices,
which is provided in a separate document (552246).



Intuitive Surgical, Inc.
1266 Kifer Road
Sunnyvale, CA 94086 USA
www.intuitivesurgical.com
Customer Service 800.876.1310

EC REP

Intuitive Surgical Sàrl
1, chemin des Mürriers
1170 Aubonne Switzerland
Customer Service (Europe)
+800.0821.2020



Copyright

© 2018 Intuitive Surgical, Inc. All rights reserved.

Product names are trademarks or registered trademarks of their respective holders.

Trademarks

Intuitive Surgical, da Vinci Xi, da Vinci X and *EndoWrist*, are trademarks or registered trademarks of Intuitive Surgical, Inc. Product names are trademarks or registered trademarks of their respective holders.

Rx Only

Table of Contents

INTRODUCTION	4
APPENDIX A	5
Validated Instruments and Accessories	
8 mm Instruments	5
BIBLIOGRAPHY	6

Introduction

This addendum is intended to be used together with the *da Vinci Xi* Reprocessing Instructions Appendices.

The contents of this addendum are supplemental instructions for Force Bipolar.

REPROCESSING INSTRUCTIONS	REPROCESSING INSTRUCTIONS PART NUMBER
<i>da Vinci Xi</i> Reprocessing Instructions Appendices	552246

Appendix A

Validated Instruments and Accessories

8 mm Instruments

da Vinci Xi PRODUCT CODES	PRODUCT DESCRIPTION	NUMBER OF USES	NUMBER OF REPROCESSING CYCLES
470405	8mm, Force Bipolar	10	15

Bibliography

A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices. AAMI TIR 30:2011.

Comprehensive guide to steam sterilization and sterility assurance in health care facilities. ANSI/AAMI ST79:2010 & A1:2010 & A2:2011 & A3:2012 & A4:2013(R)2014.

Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers. AAMI TIR 12:2010.

Guidance for Industry and Food and Drug Administration Staff. Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling. March 17, 2015

Human factors engineering for processing medical devices. AAM/TIR 55:2014.

Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices. ANSI/AAMI/ISO 17665-1:2006/(R)2013.

Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices. ANSI/AAMI ST81:2004/(R)2010.

Sterilization of medical devices. Information to be provided by the manufacturer for the processing of resterilizable medical devices. BS EN ISO 17664:2004.

Water for the reprocessing of medical devices. AAMI TIR34.



